

NAPA LP-15

Neonatal Airway Pressure Monitor

The NAPA LP-15 is a compact monitoring device that monitors mean airway pressure during CPAP/NHF therapy. It alarms when the mean airway pressure deviates from the user-selected high and low alarm limits.

- Measures pressures in both open and closed systems.
- Monitors and displays mean airway pressures from 0.0 to 35.0cm H₂O.
- Intuitive user interface provides adjustable audio and visual high and low pressure alarms.
- **High Pressure Alarm:** Provides a high pressure audible alarm when device detects pressure 0.1cm H₂O above preset high limit and the yellow LED up arrow blinks. The high alarm settings are from 6.0 to 35.0cm H₂O. Alarm resets when the desired pressure is resumed.
- **Low Pressure Alarm:** Provides a low pressure audible alarm when the device detects pressure 0.1cm H₂O below preset low limits. The yellow LED down arrow blinks when there is a low pressure alarm event. Low alarm settings are from 0.1 to 28.0cm H₂O. Alarm resets when the desired pressure is resumed.




**MONITORING FOR NEONATAL PATIENTS
ON CPAP/NHF THERAPY**

For Product Info and Demo Videos,
please visit us at www.drwmedical.com

Use the NP-CATH with your NAPAP LP-15 to MEASURE DELIVERED PRESSURES during Nasal High Flow Therapy!

- ☑ Multiple variables (flow rate, inspiratory demand, prong placement and size, etc.) create a challenging environment for delivering safe and effective NHF therapy. NAPA monitoring accounts for all these variables and measures the pressure being delivered to the patient.
- ☑ **Low flow rates result in low or negative pressures**, and can actually increase the patient's work of breathing when patient inspiratory demands are not being met.
- ☑ NHF monitoring allows clinicians to deliver **higher, evidenced-based flow rates safely**, creating a dynamic residual pressure even during flow-based therapies. High pressure alarms alert clinicians to potential over-pressure situations, thus creating an environment where therapeutically sufficient flows can be delivered safely.
- ☑ The NP-Cath accurately monitors pressures with any NHF cannula.

Measuring Pressures During NHF Therapy with the NAPA Monitor



For additional product info, supportive studies, and demo videos, please visit us at www.drwmedical.com

Distributed by:
DRW Medical

DRW Medical, LLC
2710 Concord Rd., Aston, PA 19014
Dan Tatum: 610-996-5308 Jen Myers: 610-787-2380